

Sulfur with Ultra Low Diesel Fuel, Set of 8 calibration standards, concentrations of 20, 40, 60, 80 and 100ppm, 100ml each (Standard for XRF)

Art. ID ASI-ULSD20-100-Set(5x100ML)
Unit 5 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1202 III E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-ULSD20-100-01	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	20	ppm	EPA 80.580 to 80.585 and ASTM D7212, D72 20	
ASI-ULSD20-100-02	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	40	ppm	EPA 80.580 to 80.585 and ASTM D7212, D72 20	
ASI-ULSD20-100-03	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	60	ppm	EPA 80.580 to 80.585 and ASTM D7212, D72 20	
ASI-ULSD20-100-04	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	80	ppm	EPA 80.580 to 80.585 and ASTM D7212, D72 20	
ASI-ULSD20-100-05	Sulfur (S) (as Di-n-butyl sulfide)	[544-40-1]	100	ppm	EPA 80.580 to 80.585 and ASTM D7212, D72 20	