

Lead in gasoline, 7 standards per set, 100 mL each

Art. ID ASI-PBG7C-Set(7x100ML)
Unit 7 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1203 II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|--------------------|-------------------|-------------|---------------------|-------|------------------------|--------|
| ASI-PBG7C-01-100ML | Lead (Pb) | [7439-92-1] | 0 | g/gal | ASTM D5059 , PART C | |
| ASI-PBG7C-02-100ML | Lead (Pb) | [7439-92-1] | 0,001 | g/gal | ASTM D5059 , PART C | |
| ASI-PBG7C-03-100ML | Lead (Pb) | [7439-92-1] | 0,005 | g/gal | ASTM D5059 , PART C | |
| ASI-PBG7C-04-100ML | Lead (Pb) | [7439-92-1] | 0,01 | g/gal | ASTM D5059 , PART C | |
| ASI-PBG7C-05-100ML | Lead (Pb) | [7439-92-1] | 0,05 | g/gal | ASTM D5059 , PART C | |
| ASI-PBG7C-06-100ML | Lead (Pb) | [7439-92-1] | 0,1 | g/gal | ASTM D5059 , PART C | |
| ASI-PBG7C-07-100ML | Lead (Pb) | [7439-92-1] | 0,3 | g/gal | ASTM D5059 , PART C | |