

Lead in gasoline, 7 standards per set, 100 mL each

Art. ID ASI-PBG7C-Set(7x100ML)
Unit 7 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1203 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-PBG7C-01-100ML	Lead (Pb)	[7439-92-1]	0	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-02-100ML	Lead (Pb)	[7439-92-1]	0,001	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-03-100ML	Lead (Pb)	[7439-92-1]	0,005	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-04-100ML	Lead (Pb)	[7439-92-1]	0,01	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-05-100ML	Lead (Pb)	[7439-92-1]	0,05	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-06-100ML	Lead (Pb)	[7439-92-1]	0,1	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-07-100ML	Lead (Pb)	[7439-92-1]	0,3	g/gal	ASTM D5059 , PART C	