

**Lead in gasoline. This set is made up of PBG7A and PBG7C 12 standards per set, 100 mL each**

Art. ID ASI-PBG12-Set(12x100ML)  
Unit 12 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1203 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-PBG7A-01-100ML	Lead (Pb)	[7439-92-1]	0	g/gal	ASTM D5059 , PART A	
ASI-PBG7C-02-100ML	Lead (Pb)	[7439-92-1]	0,001	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-03-100ML	Lead (Pb)	[7439-92-1]	0,005	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-04-100ML	Lead (Pb)	[7439-92-1]	0,01	g/gal	ASTM D5059 , PART C	
ASI-PBG7C-05-100ML	Lead (Pb)	[7439-92-1]	0,05	g/gal	ASTM D5059 , PART C	
ASI-PBG7A-02-100ML	Lead (Pb)	[7439-92-1]	0,1	g/gal	ASTM D5059 , PART A	
ASI-PBG7C-07-100ML	Lead (Pb)	[7439-92-1]	0,3	g/gal	ASTM D5059 , PART C	
ASI-PBG7A-03-100ML	Lead (Pb)	[7439-92-1]	1	g/gal	ASTM D5059 , PART A	
ASI-PBG7A-04-100ML	Lead (Pb)	[7439-92-1]	2	g/gal	ASTM D5059 , PART A	
ASI-PBG7A-05-100ML	Lead (Pb)	[7439-92-1]	3	g/gal	ASTM D5059 , PART A	
ASI-PBG7A-06-100ML	Lead (Pb)	[7439-92-1]	4	g/gal	ASTM D5059 , PART A	
ASI-PBG7A-07-100ML	Lead (Pb)	[7439-92-1]	5	g/gal	ASTM D5059 , PART A	