

Lead in Isooctane, 7 standards per set, 100 mL each

Art. ID ASI-PBIO7A-Set(7x100ML)
Unit 7 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-PBIO7A-01-100ML	Lead (Pb)	[7439-92-1]	0	g/gal	ASTM D5059 , PART A	
ASI-PBIO7A-02-100ML	Lead (Pb)	[7439-92-1]	0,1	g/gal	ASTM D5059 , PART A	
ASI-PBIO7A-03-100ML	Lead (Pb)	[7439-92-1]	1	g/gal	ASTM D5059 , PART A	
ASI-PBIO7A-04-100ML	Lead (Pb)	[7439-92-1]	2	g/gal	ASTM D5059 , PART A	
ASI-PBIO7A-05-100ML	Lead (Pb)	[7439-92-1]	3	g/gal	ASTM D5059 , PART A	
ASI-PBIO7A-06-100ML	Lead (Pb)	[7439-92-1]	4	g/gal	ASTM D5059 , PART A	
ASI-PBIO7A-07-100ML	Lead (Pb)	[7439-92-1]	5	g/gal	ASTM D5059 , PART A	