

**Nickel Reagent and 0.05% Manganese internal standard in xylene-mineral oil. Conc. Range: 0-0.0100wt%, 7 standards per set, 100 mL each**

Art. ID ASI-NiISO7-Set(7x100ML)  
Unit 7 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1

| Text/Information    | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method           | Source |
|---------------------|-------------------|-------------|---------------------|------|------------------|--------|
| ASI-NiISO7-01-100ML | Nickel (Ni)       | [7440-02-0] | 0                   | %    | ISO/CD 145<br>96 |        |
| ASI-NiISO7-02-100ML | Nickel (Ni)       | [7440-02-0] | 0,0005              | %    | ISO/CD 145<br>96 |        |
| ASI-NiISO7-03-100ML | Nickel (Ni)       | [7440-02-0] | 0,001               | %    | ISO/CD 145<br>96 |        |
| ASI-NiISO7-04-100ML | Nickel (Ni)       | [7440-02-0] | 0,0025              | %    | ISO/CD 145<br>96 |        |
| ASI-NiISO7-05-100ML | Nickel (Ni)       | [7440-02-0] | 0,005               | %    | ISO/CD 145<br>96 |        |
| ASI-NiISO7-06-100ML | Nickel (Ni)       | [7440-02-0] | 0,007               | %    | ISO/CD 145<br>96 |        |
| ASI-NiISO7-07-100ML | Nickel (Ni)       | [7440-02-0] | 0,01                | %    | ISO/CD 145<br>96 |        |