

Nickel Reagent and 0.05% Manganese internal standard in xylene-mineral oil. Conc. Range: 0-0.0100wt%, 7 standards per set, 100 mL each

Art. ID ASI-NiISO7-Set(7x100ML)
Unit 7 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-NiISO7-01-100ML	Nickel (Ni)	[7440-02-0]	0	%	ISO/CD 145 96	
ASI-NiISO7-02-100ML	Nickel (Ni)	[7440-02-0]	0,0005	%	ISO/CD 145 96	
ASI-NiISO7-03-100ML	Nickel (Ni)	[7440-02-0]	0,001	%	ISO/CD 145 96	
ASI-NiISO7-04-100ML	Nickel (Ni)	[7440-02-0]	0,0025	%	ISO/CD 145 96	
ASI-NiISO7-05-100ML	Nickel (Ni)	[7440-02-0]	0,005	%	ISO/CD 145 96	
ASI-NiISO7-06-100ML	Nickel (Ni)	[7440-02-0]	0,007	%	ISO/CD 145 96	
ASI-NiISO7-07-100ML	Nickel (Ni)	[7440-02-0]	0,01	%	ISO/CD 145 96	