

Vanadium reagent and 0.05% Manganese internal standard in xylene-mineral oil. Conc. Range: 0.03-0.100wt%, 7 standards per set, 100 mL each

Art. ID ASI-VISOH7-Set(7x100ML)
Unit 7 x 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-VISOH7-01-100ML	Vanadium (V)	[7440-62-2]	0	0	ISO/CD 145 96	
ASI-VISOH7-02-100ML	Vanadium (V)	[7440-62-2]	0,03	0	ISO/CD 145 96	
ASI-VISOH7-03-100ML	Vanadium (V)	[7440-62-2]	0,04	0	ISO/CD 145 96	
ASI-VISOH7-04-100ML	Vanadium (V)	[7440-62-2]	0,05	0	ISO/CD 145 96	
ASI-VISOH7-05-100ML	Vanadium (V)	[7440-62-2]	0,06	0	ISO/CD 145 96	
ASI-VISOH7-06-100ML	Vanadium (V)	[7440-62-2]	0,08	0	ISO/CD 145 96	
ASI-VISOH7-07-100ML	Vanadium (V)	[7440-62-2]	0,1	0	ISO/CD 145 96	