

**Vanadium reagent and 0.05% Manganese internal standard in xylene-mineral oil. Conc. Range: 0.0005-0.0200wt%, 9 standards per set, 100 mL each**

Art. ID ASI-VISOL9-Set(9x100ML)  
Unit 9 x 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1307 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ASI-VISOL9-01-100ML	Vanadium (V)	[7440-62-2]	0,0005	%	ISO/CD 145 96	
ASI-VISOL9-02-100ML	Vanadium (V)	[7440-62-2]	0,0025	%	ISO/CD 145 96	
ASI-VISOL9-03-100ML	Vanadium (V)	[7440-62-2]	0,005	%	ISO/CD 145 96	
ASI-VISOL9-04-100ML	Vanadium (V)	[7440-62-2]	0,0075	%	ISO/CD 145 96	
ASI-VISOL9-05-100ML	Vanadium (V)	[7440-62-2]	0,01	%	ISO/CD 145 96	
ASI-VISOL9-06-100ML	Vanadium (V)	[7440-62-2]	0,0125	%	ISO/CD 145 96	
ASI-VISOL9-07-100ML	Vanadium (V)	[7440-62-2]	0,015	%	ISO/CD 145 96	
ASI-VISOL9-08-100ML	Vanadium (V)	[7440-62-2]	0,0175	%	ISO/CD 145 96	
ASI-VISOL9-09-100ML	Vanadium (V)	[7440-62-2]	0,02	%	ISO/CD 145 96	