

**Multi-element 11 wear metals (Sulfur-Free) each 900 ppm in 20cst oil for analysis by SRF, AA, ICP or AE (for ASTM D4628, D4927, D4951, D5056, D5185, D6443, D6481)**

Art. ID ASI-WM11-900-200ML  
Unit 200 mL  
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	900	ppm		
	Chromium (Cr)	[7440-47-3]	900	ppm		
	Copper (Cu)	[7440-50-8]	900	ppm		
	Iron (Fe)	[7439-89-6]	900	ppm		
	Magnesium (Mg)	[7439-95-4]	900	ppm		
	Sodium (Na)	[7440-23-5]	900	ppm		
	Nickel (Ni)	[7440-02-0]	900	ppm		
	Lead (Pb)	[7439-92-1]	900	ppm		
	Silicon (Si)	[7440-21-3]	900	ppm		
	Tin (Sn)	[7440-31-5]	900	ppm		
	Titanium (Ti)	[7440-32-6]	900	ppm		