

Multi-element 12 wear metals (Sulfur-Free) each 10 ppm in 20cst oil for analysis by SRF, AA, ICP or AE (for ASTM D4628, D4927, D4951, D5056, D5185, D6443, D6481)

Art. ID ASI-WM12-10-200ML
Unit 200 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 10 | ppm | | |
| | Aluminium (Al) | [7429-90-5] | 10 | ppm | | |
| | Chromium (Cr) | [7440-47-3] | 10 | ppm | | |
| | Copper (Cu) | [7440-50-8] | 10 | ppm | | |
| | Iron (Fe) | [7439-89-6] | 10 | ppm | | |
| | Magnesium (Mg) | [7439-95-4] | 10 | ppm | | |
| | Sodium (Na) | [7440-23-5] | 10 | ppm | | |
| | Nickel (Ni) | [7440-02-0] | 10 | ppm | | |
| | Lead (Pb) | [7439-92-1] | 10 | ppm | | |
| | Silicon (Si) | [7440-21-3] | 10 | ppm | | |
| | Tin (Sn) | [7440-31-5] | 10 | ppm | | |
| | Titanium (Ti) | [7440-32-6] | 10 | ppm | | |