

Multi-element 20 wear metals (Sulfur-Free) each 900 ppm in 20cst oil for analysis by SRF, AA, ICP or AE (for ASTM D4628, D4927, D4951, D5056, D5185, D6443, D6481)

Art. ID ASI-WM20-900-100ML
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 900 | ppm | | |
| | Boron (B) | [7440-42-8] | 900 | ppm | | |
| | Barium (Ba) | [7440-39-3] | 900 | ppm | | |
| | Calcium (Ca) | [7440-70-2] | 900 | ppm | | |
| | Cadmium (Cd) | [7440-43-9] | 900 | ppm | | |
| | Chromium (Cr) | [7440-47-3] | 900 | ppm | | |
| | Copper (Cu) | [7440-50-8] | 900 | ppm | | |
| | Iron (Fe) | [7439-89-6] | 900 | ppm | | |
| | Magnesium (Mg) | [7439-95-4] | 900 | ppm | | |
| | Manganese (Mn) | [7439-96-5] | 900 | ppm | | |
| | Molybdenum (Mo) | [7439-98-7] | 900 | ppm | | |
| | Sodium (Na) | [7440-23-5] | 900 | ppm | | |
| | Nickel (Ni) | [7440-02-0] | 900 | ppm | | |
| | Lead (Pb) | [7439-92-1] | 900 | ppm | | |
| | Phosphorus (P) | [7723-14-0] | 900 | ppm | | |
| | Silicon (Si) | [7440-21-3] | 900 | ppm | | |
| | Tin (Sn) | [7440-31-5] | 900 | ppm | | |
| | Titanium (Ti) | [7440-32-6] | 900 | ppm | | |
| | Vanadium (V) | [7440-62-2] | 900 | ppm | | |
| | Zinc (Zn) | [7440-66-6] | 900 | ppm | | |