

**CHEMPLEX, CUPS-series 1000: TrimLess, 100 x Ø31.2mm / Height 23.6mm, Open Ends: Double**

Art. ID	CX-1060
Unit	100 units
Deliverydetails	No Dangerous Good /not restricted

**Description**

SERIES 1000: TrimLess® Sleeve Sample Cups - TrimLess® sleeve completely envelops thin - film - Eliminates thin - film trimming - No annoying clippings and annoying static cling - ThermoPlastic® seal venting - Available in Single and Double Open-Ends