

CHEMPLEX, CUPS-series 1000: TrimLess, 100 x Ø38.4mm / Height 22.9mm, Open Ends: Single

Art. ID	CX-1075
Unit	100 units
Deliverydetails	No Dangerous Good /not restricted

Description

SERIES 1000: TrimLess® Sleeve Sample Cups - TrimLess® sleeve completely envelops thin - film - Eliminates thin - film trimming - No annoying clippings and annoying static cling - ThermoPlastic® seal venting - Available in Single and Double Open-Ends