

CHEMPLEX, CUPS-series 1000: TrimLess, 100 x Ø42.9mm / Height 20.3mm, Open Ends: Double

Art. ID	CX-1080
Unit	100 units
Deliverydetails	No Dangerous Good /not restricted

Description

SERIES 1000: TrimLess® Sleeve Sample Cups - TrimLess® sleeve completely envelops thin - film - Eliminates thin - film trimming - No annoying clippings and annoying static cling - ThermoPlastic® seal venting - Available in Single and Double Open-Ends