

CHEMPLEX, CUPS-series 1300: TrimLess, 100 x Ø30.7mm / Height 22.9mm, Open Ends: Double

Art. ID	CX-1330-SE
Unit	100 units
Deliverydetails	No Dangerous Good /not restricted

Description

SERIES 1300: TrimLess® Sleeved Sample Cups - Re - sealable cap with integrated vent - External overflow reservoir integrated in cap - "Snap - On Ring" with bead - to - indent geometry - Available in nominal 32 and 40mm diameters - Serrated edged "Snap - On Rings" alternatively available to facilitate thin - film trimming