

**CHEMPLEX, CUPS-series 1800: Single Open-Ended Sample Cups, 100 x Ø31.0mm / Height 22.1mm**

|                 |                                   |
|-----------------|-----------------------------------|
| Art. ID         | CX-1830-SE                        |
| Unit            | 100 units                         |
| Deliverydetails | No Dangerous Good /not restricted |

**Description**

SERIES 1800: Single Open - Ended Sample Cups - External overflow reservoir - Snap - On Ring thin - film attachment - ThermoPlastic® seal venting