

CHEMPLEX, CUPS-series 1850: SpectroSulfur, 100 x Ø42.9mm / Height 19.6mm, Open Ends: Single - please ask for availability

Art. ID	CX-1850
Unit	100 units
Deliverydetails	No Dangerous Good /not restricted

Description

1850: SpectroSulfur® Analyzer Low Profile XRF Sample Cup - Single Open End - Shallow single open - ended - External overflow reservoir - "Snap - On Ring" thin - film attachment - ThermoPlastic® seal venting - For XOS Sindie® Systems and other Sulfur Analyzers