

**CHEMPLEX, CUPS-series 1935-OX: XRF Sample Cups and Vented Cap Replacements for Oxford Analyzers, Double Open Ends, 100 x Ø31.2mm / Height 38.4mm**

Art. ID	CX-1935-OX
Unit	100 units
Deliverydetails	No Dangerous Good /not restricted

Description

1935 - OX: XRF Sample Cups and Vented Cap Replacements for Oxford Analyzers - Double Open Ends - Replacement for Oxford 54 - CK - 10 - Includes vented friction fitting cover - External overflow reservoir