

**L-Arginine,PharmaGrade, Ajinomoto, EP, USP, manufactured under appropriate GMP controls for Pharma or Biopharmaceutical production, suitable for cell culture (SIGMA)**

Art. ID                      SAF-A4474-1KG  
Unit                              1 kg  
Deliverydetails              No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | L-Arginine        | [74-79-3]  |                     |      |        |        |