

Copper base, E - Cu, Setting-up-sample, disc Ø 40 mm x 25 mm

Art. ID SUS-2.0070
Unit disc
Deliverydetails No Dangerous Good /not restricted

Description

Please note that the values varies from lot to lot.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Zinc (Zn)	[7440-66-6]	<0,001	%		
	Lead (Pb)	[7439-92-1]	<0,001	%		
	Tin (Sn)	[7440-31-5]	<0,001	%		
	Phosphorus (P)	[7723-14-0]	<0,001	%		
	Manganese (Mn)	[7439-96-5]	<0,001	%		
	Iron (Fe)	[7439-89-6]	<0,001	%		
	Nickel (Ni)	[7440-02-0]	<0,001	%		
	Silicon (Si)	[7440-21-3]	<0,001	%		
	Chromium (Cr)	[7440-47-3]	<0,001	%		
	Antimony (Sb)	[7440-36-0]	<0,001	%		
	Aluminium (Al)	[7429-90-5]	<0,001	%		