

## **Certified Reference Material - Sulfide, waste water (WasteWatR TM)**

Art. ID                      ERA-071  
Unit                         each  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1

### Description

One 10 mL flame-sealed ampule yields up to 1 liter after dilution. Preserved sample is guaranteed stable. Analyze for sulfide by titrimetric or colorimetric methods or ISE. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfide		2 - 10	mg/L		