

Certified Reference Material - PAHs on Polyurethane foam

Art. ID	ERA-1113
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

One 2 mL flame-sealed ampule plus one polyurethane foam (PUF). Contains at least 13 analytes, randomly selected from the list below, at 10?200 µg/sample after preparation. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene	[83-32-9]	10 - 200	µg/sample		
	Acenaphthylene	[208-96-8]	10 - 200	µg/sample		
	Anthracene	[120-12-7]	10 - 200	µg/sample		
	Benzo(a)anthracene	[56-55-3]	10 - 200	µg/sample		
	Benzo(a)pyrene	[50-32-8]	10 - 200	µg/sample		
	Benzo(b)fluoranthene	[205-99-2]	10 - 200	µg/sample		
	Benzo(g,h,i)perylene	[191-24-2]	10 - 200	µg/sample		
	Benzo(k)fluoranthene	[207-08-9]	10 - 200	µg/sample		
	Chrysene	[218-01-9]	10 - 200	µg/sample		
	Dibenz(a,h)anthracene	[53-70-3]	10 - 200	µg/sample		
	Fluoranthene	[206-44-0]	10 - 200	µg/sample		
	Fluorene	[86-73-7]	10 - 200	µg/sample		
	Indeno(1,2,3-cd)pyrene	[193-39-5]	10 - 200	µg/sample		
	Naphthalene	[91-20-3]	10 - 200	µg/sample		
	Phenanthrene	[85-01-8]	10 - 200	µg/sample		
	Pyrene	[129-00-0]	10 - 200	µg/sample		