

Certified Reference Material - Metals on Filter Paper

Art. ID	ERA-1125
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

One filter paper sample packaged in a 50 mm polystyrene petri dish containing a single 47 mm tissue quartz filter ready for use with EPA Method 29 or other applicable methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	25 - 250	µg/filter		
	Arsenic (As)	[7440-38-2]	20 - 250	µg/filter		
	Barium (Ba)	[7440-39-3]	20 - 250	µg/filter		
	Beryllium (Be)	[7440-41-7]	10 - 250	µg/filter		
	Cadmium (Cd)	[7440-43-9]	10 - 250	µg/filter		
	Chromium (Cr)	[7440-47-3]	15 - 250	µg/filter		
	Cobalt (Co)	[7440-48-4]	10 - 250	µg/filter		
	Copper (Cu)	[7440-50-8]	10 - 250	µg/filter		
	Lead (Pb)	[7439-92-1]	20 - 350	µg/filter		
	Manganese (Mn)	[7439-96-5]	10 - 250	µg/filter		
	Nickel (Ni)	[7440-02-0]	20 - 250	µg/filter		
	Phosphorous		10 - 250	µg/filter		
	Selenium (Se)	[7782-49-2]	20 - 250	µg/filter		
	Silver (Ag)	[7440-22-4]	30 - 250	µg/filter		
	Thallium (Tl)	[7440-28-0]	30 - 250	µg/filter		
	Zinc (Zn)	[7440-66-6]	20 - 250	µg/filter		