

## **Certified Reference Material - Sulfur Dioxide in Impinger Solution**

Art. ID                      ERA-1143  
Unit                         each  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1830 II E2

### Description

One impinger solution sample packaged in a 15 mL screw-top vial containing approximately 14 mL of standard concentrate for use with EPA Method 6 and Method 8, or other applicable methods. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Sulfur Dioxide	[7446-09-5]	50 - 2000	mg/dscm		