

## **Certified Reference Material - Ammonia in Impinger Solution**

Art. ID                      ERA-1145  
Unit                         each  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

### Description

One impinger solution sample packaged in a 15 mL screw-top vial containing approximately 14 mL of standard concentrate for use with EPA CTM 027, or other applicable methods. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonium (NH <sub>4</sub> <sup>+</sup> )	[14798-03-9]	0,1 - 10	mg/L		