

## **Certified Reference Material - Solids, Concentrate, waste water (WasteWatR TM)**

Art. ID	ERA-4032
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 23 mL screw-cap vial with a powder yields 1 liter after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		20 - 100	mg/L		
	Total dissolved solids at 180 °C		140 - 800	mg/L		
	Total solids at 105 °C		140 - 800	mg/L		