

Certified Reference Material - Oil & Grease Concentrate, waste water (WasteWatR TM)

Art. ID	ERA-4122
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1274 II E2

Description

One 24 mL screw-cap vial yields up to 2 liters after dilution. Gravimetric analysis only. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Hexane Extractable Material (O and G) (Gravimetric)		20 - 200	mg/L		