

Certified Reference Material - Solids, waste water (WasteWatR TM)

Art. ID	ERA-499
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

One 500 mL whole-volume bottle is ready to analyze. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total suspended solids (TSS)		20 - 100	mg/L		
	Total dissolved solids at 180 °C		140 - 800	mg/L		
	Total solids at 105 °C		140 - 800	mg/L		