

## **Certified Reference Material - Total Residual Chlorine, waste water (WasteWatR TM)**

Art. ID	ERA-501
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Use with titrimetric or colorimetric methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total residual chlorine		0,5 - 3	mg/L		
	Free Residual Chlorine		0,5 - 3	mg/L		