

Certified Reference Material - Cyanide, waste water (WasteWatR TM)

Art. ID	ERA-502
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1

Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total cyanide		0,1 - 1	mg/L		
	Amenable Cyanide		0,1 - 1	mg/L		
	Available Cyanide		0,1 - 1	mg/L		