

## **Certified Reference Material - Cyanide, waste water (WasteWatR TM)**

|                 |   |
|-----------------|---|
| Art. ID         | ERA-502   |
| Unit            | each  |
| Deliverydetails | ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1 |

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | Total cyanide     |            | 0,1 - 1             | mg/L |        |        |
|                  | Amenable Cyanide  |            | 0,1 - 1             | mg/L |        |        |
|                  | Available Cyanide |            | 0,1 - 1             | mg/L |        |        |