

## **Certified Reference Material - Simple Nutrients, waste water (WasteWatR TM)**

Art. ID ERA-505  
Unit each  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonia (as N)		1,0 - 20,0	mg/L		
	Nitrate + Nitrite as N		2,5 - 25	mg/L		
	Nitrate (NO <sub>3</sub> -) as N		2,0 - 25,0	mg/L		
	Ortho-Phosphate as P		0,5 - 5,5	mg/L		
	Total nitrogen		3,0 - 45,0	mg/L		