

Certified Reference Material - Mercury, waste water (WasteWatR TM)

Art. ID ERA-514
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

One 15 mL screw-cap vial yields up to 1 liter after dilution. Analyze for total mercury. The certified values vary from batch to batch. Please ask for the certificate.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Mercury (Hg) | [7439-97-6] | 3,0 - 30,0 | µg/L | | |