

## **Certified Reference Material - Solids, Concentrate, potable water (PotableWatR TM)**

Art. ID	ERA-5152
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 23 mL screw-cap vial with a powder yields 1 liter after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Non-filterable residue (TSS)		23 - 100	mg/L		
	Total dissolved solids at 180 °C		100 - 1000	mg/L		
	Total solids at 105 °C		123 - 1100	mg/L		