

Certified Reference Material - Demand, waste water (WasteWatR TM)

Art. ID	ERA-516
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Biochemical Oxygen Demand (BOD)		18 - 230	mg/L		
	Carbonous BOD (CBOD)		18 - 230	mg/L		
	Chemical Oxygen Demand (COD)		30 - 250	mg/L		
	Total Organic Carbon (TOC)		6 - 100	mg/L		