

## **Certified Reference Material - Tin & Titanium, waste water (WasteWatR TM)**

Art. ID                      ERA-517  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

One 15 mL screw-cap vial yields up to 1 liter after dilution. Use with AA, ICP-OES or ICP-MS methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tin (Sn)	[7440-31-5]	1000 - 5000	µg/L		
	Titanium (Ti)	[7440-32-6]	80 - 300	µg/L		