

Certified Reference Material - 1 liter Oil & Grease, waste water (WasteWatR TM)

Art. ID	ERA-518
Unit	1 L
Deliverydetails	No Dangerous Good /not restricted

Description

One liter whole-volume glass bottle with a 33-430 thread is ready to analyze. For SPE compatible QC (33-400 thread), use Cat #818 1 liter Boston Round Oil & Grease. Please specify if for SPE compatible PT. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	n-Hexane Extractable Material(O&G)(Grav)		20 - 200	mg/L		
	n-Hexane Extractable Material(O&G)(IR)		20 - 200	mg/L		