

Certified Reference Material - HEM / SGT-HEM, waste water (WasteWatR TM)

Art. ID ERA-519
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

Description

One 5 mL flame-sealed ampule yields up to 2 liters after dilution. Contains both hexadecane and stearic acid. Note: If a NELAC compliant PT is required, use Cat. #582 or Cat. #4120. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	HEM		5 - 100	mg/L		
	SGT-HEM		5 - 100	mg/L		