

## **Certified Reference Material - Complex Nutrients, waste water (WasteWatR TM)**

Art. ID	ERA-525
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

### Description

One 15 mL screw-cap vial yields up to 2 liters after dilution. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Kjeldahl-Nitrogen as N		3,0 - 35,0	mg/L		
	Total Phosphorus as P		0,5 - 10	mg/L		