

Certified Reference Material - Inorganic Disinfection 1 QC

Art. ID	ERA-5272
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

One 23 mL screw-cap vial yields up to 4 liters after dilution. The diluted standard will contain chlorate at 60 to 180 µg/L. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chlorate		60 - 180	µg/L		
	Chlorite		100 - 1000	µg/L		