

Certified Reference Material - Mercury, potable water (PotableWatR TM)

Art. ID ERA-666
Unit each
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

One 15 mL screw-cap vial yields up to 1 liter after dilution. Use with CVAA, ICP-MS or CVAFS methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mercury (Hg)	[7439-97-6]	0,5 - 10	µg/L		