

## **Certified Reference Material - Inorganics, potable water (PotableWatR TM)**

Art. ID	ERA-698
Unit	each
Deliverydetails	No Dangerous Good /not restricted

### Description

One 500 mL whole-volume bottle is ready to analyze. The CRM is also certified for Sodium. For a Sodium PT, order Hardness, Cat. #555. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Alkalinity as CaCO <sub>3</sub>		25 - 200	mg/L		
	Chloride (Cl)		20 - 160	mg/L		
	Fluoride (F)		1 - 8	mg/L		
	Nitrate + Nitrite as N		3 - 10	mg/L		
	Nitrate (NO <sub>3</sub> -) as N		3 - 10	mg/L		
	Potassium (K)	[7440-09-7]	10 - 40	mg/L		
	Sulfate		25 - 250	mg/L		
	Specific Conductance at 25°C		130 - 1300	µmhos/cm		
	Total filterable residu e (TDS) at 180°C		100- 1000	mg/L		