

Certified Reference Material - Toxaphene, potable water (PotableWatR TM)

Art. ID ERA-700
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains toxaphene at 2-20 µg/L. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Campechlor (Toxaphene) (total)		3,0 -20,0	µg/L		