

### **Certified Reference Material - Chlordane, potable water (PotableWatR TM)**

Art. ID                      ERA-705  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

#### Description

One 2 mL flame-sealed ampule yields up to 2 liters after dilution. Each standard contains technical chlordane at 2-20 µg/L. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Chlordane (total)		2,0 - 20,0	µg/L		