

Certified Reference Material - Toxaphene, priority pollutant (PriorityPollutnT TM)

Art. ID ERA-717
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One 2 mL flame-sealed ampule yields up to 2 liters of sample after dilution. Contains toxaphene at 20-100 µg/L. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Campechlor (Toxaphene)	[8001-35-2]	20 - 100	µg/L		