

Certified Reference Material - TCLP Semivolatiles

Art. ID	ERA-737
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

One 2 mL flame-sealed ampule containing a subset of the analytes listed below, each at a concentration of 0.5-2.0 mg/L after dilution. All unspiked analytes are certified at <0.5 mg/L. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,5-Trichlorophenol	[95-95-4]	0,1 - 2,0	mg/L		
	2,4,6-Trichlorophenol	[88-06-2]	0,1 - 2,0	mg/L		
	2-Methylphenol (o-Cresol)	[95-48-7]	0,1 - 2,0	mg/L		
	3- and 4-Methylphenol		0,1 - 2,0	mg/L		
	Pentachlorophenol	[87-86-5]	0,1 - 2,0	mg/L		