

### **Certified Reference Material - Silica, waste water (WasteWatR TM)**

Art. ID                      ERA-775  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

#### Description

One 60 mL poly bottle yields up to 1 liter after dilution. Analyze for silica as SiO<sub>2</sub> with colorimetric or ICP methods. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silica as SiO <sub>2</sub>		50 - 250	mg/L		