

Certified Reference Material - Turbidity, waste water (WasteWatR TM)

Art. ID ERA-777
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One 24 mL screw-cap vial yields up to 1 liter after dilution. Use with nephelometric methods. The certified values vary from batch to batch. Please ask for the certificate.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| | Turbidity | | 2,0 - 30,0 | NTU | | |