

Certified Reference Material - PCBs in Water - 1016

Art. ID ERA-860
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

One flame-sealed ampule with 2 mL of concentrate makes 1 liter after dilution. To be analyzed for Aroclor 1016 at 1-15 µg/L. The certified values vary from batch to batch. Please ask for the certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aroclor 1016	[12674-11-2]	1,0 - 15	µg/L		